

Anti-Human/Mouse TNC Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08313

3th Edition

Description: This gene encodes an extracellular matrix protein with a spatially and temporally restricted tissue distribution. This protein is homo-hexameric with disulfide-linked subunits, and contains multiple EGF-like and fibronectin type-III domains. It is implicated in guidance of migrating neurons as well as axons during development, synaptic plasticity, and neuronal regeneration.

Antigen: Synthetic peptide of human TNC

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: