

**Anti-Human/Mouse RAB24 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08359**

---

3th Edition

**Description:**RAB24 is a small GTPase of the Rab subfamily of Ras-related proteins that regulate intracellular protein trafficking. May be involved in autophagy-related processes. Unlike other Rab family members, does not interact with GDP dissociation inhibitors (GDIs), including ARHGDI A and ARHGDI B. Interacts with ZFYVE20.

**Antigen:**Synthetic peptide of human RAB24

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**