

**Anti-Human/Mouse ARHGAP5 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08369**

---

3th Edition

**Description:** Rho GTPase activating protein 5 negatively regulates RHO GTPases, a family which Human, Mouse mediate cytoskeleton changes by stimulating the hydrolysis of bound GTP. Two transcript variants encoding different isoforms have been found for this gene. GTPase-activating protein for Rho family members. May play a role in the reduction of the p21rasGTPase-activating potential of RASA1/p120GAP.

**Antigen:** Recombinant protein of human ARHGAP5

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in kidney, brain, liver and lung.