

Anti-Human ABI3BP Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08396

3th Edition

Description:TARSH gene expression is dramatically induced in mouse embryonic fibroblasts (MEFs) replicative senescence and suppressed in human lung carcinoma specimens and thyroid carcinomas. The experiment observed that ABI3 expression is reduced or lost in most carcinomas but also that there is a positive correlation between ABI3 and ABI3BP expression. And it express high levels in brain, heart, lung, liver, pancreas kidney and placenta.

Antigen:Synthetic peptide of human ABI3BP

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in brain, heart, lung, liver, pancreas kidney and placenta.