

Anti-Human/Mouse/Rat ADHFE1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08417

3th Edition

Description:The ADHFE1 gene encodes hydroxyacid-oxoacid transhydrogenase (EC 1.1.99.24), which is responsible for the oxidation of 4-hydroxybutyrate in mammalian tissues.

Antigen:Synthetic peptide of human ADHFE1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Only expressed in adult liver.1 Publication