

**Anti-Human/Mouse PLXNA1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08439**

---

3th Edition

**Description:** Plexins are a family of large, transmembrane receptors for multiple classes of semaphorins in vertebrates. Plexins are widely expressed, and regions of their extracellular domain are homologous to both scatter factor receptors and semaphorin domains.

**Antigen:** Synthetic peptide of human PLXNA1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Detected in fetal brain, lung, liver and kidney. 1 Publication