

Anti-Human/Mouse PRDM14 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08447

3th Edition

Description: This gene encodes a member of the PRDI-BF1 and RIZ homology domain containing (PRDM) family of transcriptional regulators. The encoded protein may possess histone methyltransferase activity and plays a critical role in cell pluripotency by suppressing the expression of differentiation marker genes. Expression of this gene may play a role in breast cancer

Antigen: Recombinant protein of human PRDM14

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in embryonic stem cells. Tends to be overexpressed in breast cancer (at protein level).