

**Anti-Human/Mouse ALYREF Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08460**

---

3th Edition

**Description:** The protein encoded by this gene is a heat stable, nuclear protein and functions as a molecular chaperone. It is thought to regulate dimerization, DNA binding, and transcriptional activity of basic region-leucine zipper (bZIP) proteins.

**Antigen:** Synthetic peptide of human ALYREF

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in a wide variety of cancer types.