

Anti-Human PARP10 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08541

3th Edition

Description: Poly(ADP-ribose) polymerases (PARPs), such as PARP10, regulate gene transcription by altering chromatin organization by adding ADP-ribose to histones. PARPs can also function as transcriptional cofactors.

Antigen: Synthetic peptide of human PARP10

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highly expressed in spleen and thymus. Intermediate levels in liver, kidney, pancreas, prostate, testis, ovary, intestine, and leukocytes. Low expression in heart, brain, placenta, lung, skeletal muscle, and colon.¹
Publication