

**Anti-Human/Mouse/Rat SGK2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08555**

---

3th Edition

**Description:** This gene encodes a serine/threonine protein kinase. Although this gene product is similar to serum- and glucocorticoid-induced protein kinase (SGK), this gene is not induced by serum or glucocorticoids. This gene is induced in response to signals that activate phosphatidylinositol 3-kinase, which is also true for SGK. Alternative splicing results in multiple transcript variants

**Antigen:** Synthetic peptide of human SGK2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Highly expressed in liver, kidney and pancreas, and at lower levels in brain.