

Anti-Human/Mouse/Rat TPR Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08613

3th Edition

Description: This gene encodes a large coiled-coil protein that forms intranuclear filaments attached to the inner surface of nuclear pore complexes (NPCs). The protein directly interacts with several components of the NPC. It is required for the nuclear export of mRNAs and some proteins. Oncogenic fusions of the 5' end of this gene with several different kinase genes occur in some neoplasias.

Antigen: Synthetic peptide of human TPR

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highest in testis, lung, thymus, spleen and brain, lower levels in heart, liver and kidney.