

Anti-Human/Mouse DCAF6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08619

3th Edition

Description:NRIP (Nuclear receptor interaction protein), also known as IQWD1 (IQ motif and WD repeat-containing protein 1), MSTP055, ARCAP or PC326, is an 860 amino acid protein that localizes to the nucleus and contains one IQ domain and seven WD-repeats. Expressed in testis, skeletal muscle, prostate and heart, NRIP functions as a ligand-dependent coactivator of nuclear receptors and specifically enhances the transcriptional activity of AR (androgen receptor) and GR (glucocorticoid receptor). NRIP exists as three isoforms that are produced by alternative splicing events.

Antigen:Synthetic peptide of human DCAF6

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Highly expressed in skeletal muscle and testis. Expressed to a lesser degree in heart, prostate, and adrenal gland.