

Anti-Human PNOC Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08650

3th Edition

Description: Nociception, a pain response mechanism, occurs in response to stimuli that threaten the integrity of an organism. The first synapses produced as a result of the initiation of nociception are modulated by excitatory amino acids (glutamate and aspartate) and many peptides (substance P, CGRP, CCK, endogenous opioids). Nociceptin (also designated orphanin FQ) is a neuronal peptide that is similar to opioid peptides. Nociceptin activates KOR-3 (kappa-type opioid receptor, also designated ORL1), a G protein-coupled receptor. Although similar to dynorphin A, a kappa opioid peptide, nociceptin functions to make animals hyperreactive to nociceptive stimulations. Nociceptin is also involved in locomotor behavior and may be involved in the modulation of synaptic plasticity in learning and memory.

Antigen: Synthetic peptide of human PNOC

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Predominantly expressed in the brain and spinal cord. Also expressed and secreted by peripheral blood neutrophils following degranulation.