

Anti-Human NLRP9 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08658

3th Edition

Description:The protein encoded by this gene belongs to the NALP protein family. Members of the NALP protein family typically contain a NACHT domain, a NACHT-associated domain (NAD), a C-terminal leucine-rich repeat (LRR) region, and an N-terminal pyrin domain (PYD). This protein may play a regulatory role in the innate immune system as similar family members belong to the signal-induced multiprotein complex, the inflammasome, that activates the pro-inflammatory caspases, caspase-1 and caspase-5.

Antigen:Synthetic peptide of human NLRP9

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in ileum intestinal epithelial cells. Not detected in peripheral blood mononuclear cells (PubMed:28636595). Expressed in cerebral endothelial cells and, at much lower levels, in brain pericytes (PubMed:28432035).