

Anti-Human CLEC16A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08663

3th Edition

Description: This gene encodes a member of the C-type lectin domain containing family. Single nucleotide polymorphisms in introns of this gene have been associated with diabetes mellitus, multiple sclerosis and rheumatoid arthritis. Multiple transcript variants encoding different isoforms have been found for this gene.

Antigen: Synthetic peptide of human CLEC16A

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Almost exclusively expressed in immune cells, including dendritic cells, B-lymphocytes and natural killer cells.1 Publication