

**Anti-Human/Mouse RAB39B Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08671**

---

3th Edition

**Description:** This gene encodes a member of the Rab family of proteins. Rab proteins are small GTPases that are involved in vesicular trafficking. Mutations in this gene are associated with X-linked mental retardation.

**Antigen:** Synthetic peptide of human RAB39B

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Highly expressed in the brain.1 Publication