

Anti-Human/Mouse/Rat NDRG3 Polyclonal Antibody
Polyclonal Antibody

Cat.NO.: PA08679

3th Edition

Description:NDRG3 belongs to the N-myc down-regulated gene (NDRG) family that contains 4 paralogs: NDRG1, -2, -3 and -4. The function of NDRG3 and its relationship to cancer has not been studied. NDRG3 expression is enriched in testis and prostate using gene expression data derived from massively parallel signature sequencing from 33 different human organs. NDRG3 is also expressed in both epithelial prostate cancer cells and prostatic stromal cells at both mRNA and proteins levels.

Antigen:Synthetic peptide of human NDRG3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous. Highly expressed in brain.