

Anti-Human NEUROG1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08723

3th Edition

Description: Neurogenin-1 is a protein that in humans is encoded by the NEUROG1 gene. Acts as a transcriptional regulator. Involved in the initiation of neuronal differentiation. Activates transcription by binding to the E box (5'-CANNTG-3'). Associates with chromatin to enhancer regulatory elements in genes encoding key transcriptional regulators of neurogenesis. NEUROG1 has been shown to interact with CREB-binding protein and decapentaplegic homolog 1.

Antigen: Synthetic peptide of human NEUROG1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expression restricted to the embryonic nervous system.