

Anti-Human/Mouse/Rat ADCY3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08771

3th Edition

Description: This gene encodes adenylyl cyclase 3 which is a membrane-associated enzyme and catalyzes the formation of the secondary messenger cyclic adenosine monophosphate (cAMP). This protein appears to be widely expressed in various human tissues and may be involved in a number of physiological and pathophysiological metabolic processes.

Antigen: Recombinant protein of human ADCY3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in brain, heart, kidney, liver, lung, pancreas, placenta, and skeletal muscle.