

**Anti-Human/Mouse ORM2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08774**

---

3th Edition

**Description:** This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression.

**Antigen:** Recombinant protein of human ORM2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed by the liver and secreted in plasma.