

**Anti-Human MAP1LC3C Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08871**

---

3th Edition

**Description:**Autophagy is a highly regulated bulk degradation process that plays an important role in cellular maintenance and development. MAP1LC3C is an ortholog of the yeast autophagosome protein Atg8

**Antigen:**Synthetic peptide of human MAP1LC3C

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Most abundant in placenta, lung and ovary.