

Anti-Human/Mouse/Rat SV2A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08945

3th Edition

Description: This protein plays a role in the control of regulated secretion in neural and endocrine cells, enhancing selectively low-frequency neurotransmission. Positively regulates vesicle fusion by maintaining the readily releasable pool of secretory vesicles. It is identified as the brain binding-site for the antiepileptic drug levetiracetam/lev.

Antigen: Synthetic peptide of human SV2A

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: