

Anti-Human/Mouse/Rat CAPZA3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08950

3th Edition

Description: This gene encodes an actin capping protein, one of the F-actin capping protein alpha subunit family. The encoded protein is predominantly localized to the neck region of ejaculated sperm, other immunohistochemical signals were found in the tail and postacrosomal regions. The encoded protein may also form heterodimers of alpha and beta subunits. This protein may be important in determining sperm architecture and male fertility.

Antigen: Recombinant protein of human CAPZA3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed exclusively in testis and sperm. Highest expression is found in the neck region of ejaculated sperm with lower levels found in the tail and postacrosome region.