

Anti-Human/Mouse/Rat CDC37L1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09022

3th Edition

Description: CDC37L1 is a cytoplasmic phosphoprotein that exists in complex with HSP90 (HSPCA; MIM 140571) as well as several other proteins involved in HSP90-mediated protein folding

Antigen: Recombinant protein of human CDC37L1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in brain, heart, kidney, liver, placenta and skeletal muscle.