

Anti-Human IFNA16 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09034

3th Edition

Description: Interferon, alpha 16 is produced by macrophages and has antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

Antigen: Synthetic peptide of human IFNA16

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: