

Anti-Human/Mouse HRASLS2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09075

3th Edition

Description:HRASLS2 (HRAS-like suppressor 2), also known as PLA1/2-2, is a 162 amino acid Cytoplasmic protein that belongs to the H-rev107 family. Expressed in liver, kidney, small intestine and colon, HRASLS2 exhibits PLA1/2 activity by catalyzing the calcium-independent hydrolysis of acyl groups in various phosphatidylcholines (PCs) and phosphatidylethanolamine (PE).

Antigen:Synthetic peptide of human HRASLS2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in liver, kidney, small intestine testis and colon (PubMed:19615464). Undetectable in testis, placenta, salivary gland and fetal brain (PubMed:18163183).