

Anti-Human/Mouse GREB1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09121

3th Edition

Description: This gene is an estrogen-responsive gene that is an early response gene in the estrogen receptor-regulated pathway. It is thought to play an important role in hormone-responsive tissues and cancer. Three alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

Antigen: Synthetic peptide of human GREB1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in proliferating prostatic tissue and prostate cancer.