

Anti-Human ADCY1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09131

3th Edition

Description: This gene encodes a form of adenylate cyclase expressed in brain. A similar protein in mouse is involved in pattern formation of the brain.

Antigen: Synthetic peptide of human ADCY1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Brain, retina and adrenal medulla.