

**Anti-Human/Mouse/Rat ADHFE1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09138**

---

3th Edition

**Description:**The ADHFE1 gene encodes hydroxyacid-oxoacid transhydrogenase (EC 1.1.99.24), which is responsible for the oxidation of 4-hydroxybutyrate in mammalian tissues (Kardon et al., 2006 [PubMed 16616524]).

**Antigen:**Synthetic peptide of human ADHFE1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Only expressed in adult liver.1 Publication