

Anti-Human/Mouse/Rat Glucagon Receptor Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09158

3th Edition

Description: The protein encoded by this gene is a glucagon receptor that is important in controlling blood glucose levels. Defects in this gene are a cause of non-insulin-dependent diabetes mellitus (NIDDM).

Antigen: Synthetic peptide of human GCGR

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: