

**Anti-Human/Mouse APOBEC3D Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09193**

---

3th Edition

**Description:** This gene is a member of the cytidine deaminase gene family. It is one of a group of related genes found in a cluster, thought to result from gene duplication, on chromosome 22. Members of the cluster encode proteins that are structurally and functionally related to the C to U RNA-editing cytidine deaminase APOBEC1 and inhibit retroviruses, such as HIV, by deaminating cytosine residues in nascent retroviral cDNA.

**Antigen:** Synthetic peptide of human APOBEC3D

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in lymphoid organs. Also detected in non-lymphoid tissues including lung.1 Publication