

Anti-Human FMN2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09220

3th Edition

Description: Formin homology (FH) domain proteins (see FMN1; MIM 136535) play a role in cytoskeletal organization and/or establishment of cell polarity.

Antigen: Synthetic peptide of human FMN2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed almost exclusively in the developing and mature central nervous system.