

Anti-Human/Mouse AMER1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09269

3th Edition

Description:The protein encoded by this gene upregulates transcriptional activation by the Wilms tumor protein and interacts with Human, Mouse any other proteins, including CTNNB1, APC, AXIN1, and AXIN2. Defects in this gene are a cause of osteopathia striata with cranial sclerosis (OSCS).

Antigen:Synthetic peptide of human AMER1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Detected in fetal and adult kidney, brain and spleen.1 Publication