

Anti-Human/Mouse DHRS3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09289

3th Edition

Description: Short-chain dehydrogenases/reductases (SDRs), such as DHRS3, catalyze the oxidation/reduction of a wide range of substrates, including retinoids and steroids.

Antigen: Recombinant protein of human DHRS3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed with highest levels found in heart, placenta, lung, liver, kidney, pancreas, thyroid, testis, stomach, trachea and spinal cord. Lower levels found in skeletal muscle, intestine and lymph node. No expression detected in brain. In the retina, expressed in cone but not rod outer segments.