

**Anti-Human/Mouse/Rat EPO Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09309**

---

3th Edition

**Description:** This gene is a member of the EPO/TPO family and encodes a secreted, glycosylated cytokine composed of four alpha helical bundles. The protein is found in the plasma and regulates red cell production by promoting erythroid differentiation and initiating hemoglobin synthesis. This protein also has neuroprotective activity against a variety of potential brain injuries and antiapoptotic functions in several tissue types.

**Antigen:** Synthetic peptide of human EPO

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Produced by kidney or liver of adult mammals and by liver of fetal or neonatal mammals.