

Anti-Human/Mouse DUSP23 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09320

3th Edition

Description:Dual specificity protein phosphatase 23, also known as low molecular mass dual specificity phosphatase 3 (LDP-3), is an enzyme that in humans is encoded by the DUSP23 gene.

Antigen:Recombinant protein of human DUSP23

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Widely expressed. Highly expressed in spleen, prostate, colon, adrenal gland, mammary gland, thyroid and trachea. Expressed at lower level in uterus, small intestine, bladder, bone marrow, brain, spinal cord and stomach.