

Anti-Human NPAP1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09328

3th Edition

Description: This gene is located in the Prader-Willi syndrome region on chromosome 15. This gene is biallelically expressed in adult testis and brain but is paternally imprinted in fetal brain. Defects in this gene may be associated with Prader-Willi syndrome.

Antigen: Synthetic peptide of human NPAP1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Testis-specific in adults. In fetal brain expressed only from the paternal allele.