

Anti-Human/Mouse/Rat MAPRE3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09329

3th Edition

Description:The protein encoded by this gene is a member of the RP/EB family of genes. The protein localizes to the cytoplasmic microtubule network and binds APCL, a homolog of the adenomatous polyposis coli tumor suppressor gene.

Antigen:Recombinant protein of human MAPRE3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Predominantly expressed in brain and muscle.