

**Anti-Human/Mouse EDIL3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09347**

---

3th Edition

**Description:**The protein encoded by this gene is an integrin ligand. It plays an important role in mediating angiogenesis and may be important in vessel wall remodeling and development. It also influences endothelial cell behavior.

**Antigen:**Recombinant protein of human EDIL3

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**