

**Anti-Human STK17A Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09351**

---

3th Edition

**Description:** This gene is a member of the DAP kinase-related apoptosis-inducing protein kinase family and encodes an autophosphorylated nuclear protein with a protein kinase domain. The protein has apoptosis-inducing activity.

**Antigen:** Synthetic peptide of human STK17A

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Highly expressed in placenta. Lower levels in heart, lung, skeletal muscle, kidney and pancreas.