

Anti-Human/Mouse EPHA10 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09389

3th Edition

Description: Ephrin receptors, the largest subfamily of receptor tyrosine kinases (RTKs), and their ephrin ligands are important mediators of cell-cell communication regulating cell attachment, shape, and mobility in neuronal and epithelial cells

Antigen: Recombinant protein of human EPHA10

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Mainly expressed in testis.