

**Anti-Human/Mouse/Rat MYL12B Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09467**

---

3th Edition

**Description:**The activity of nonmuscle myosin II (see MYH9; MIM 160775) is regulated by phosphorylation of a regulatory light chain, such as MRLC2. This phosphorylation results in higher MgATPase activity and the assembly of myosin II filaments

**Antigen:**Recombinant protein of human MYL12B

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Ubiquitously expressed in various hematopoietic cells.