

**Anti-Human CD80 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09481**

---

3th Edition

**Description:**The protein encoded by this gene is a membrane receptor that is activated by the binding of CD28 or CTLA-4. The activated protein induces T-cell proliferation and cytokine production. This protein can act as a receptor for adenovirus subgroup B and may play a role in lupus neuropathy.

**Antigen:**Synthetic peptide of human CD80

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed on activated B-cells, macrophages and dendritic cells.