

Anti-Human/Mouse ADGRB2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09487

3th Edition

Description:BAI1, a p53-target gene, encodes brain-specific angiogenesis inhibitor, a seven-span transmembrane protein and is thought to be a member of the secretin receptor family. Brain-specific angiogenesis proteins BAI2 and BAI3 are similar to BAI1 in structure, have similar tissue specificities and may also play a role in angiogenesis.

Antigen:Synthetic peptide of human BAI2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Strongly expressed in brain. Also detected in heart, thymus, skeletal muscle, and different cell lines.