

**Anti-Human TSPYL2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09495**

---

3th Edition

**Description:** This gene encodes a member of the testis-specific protein Y-encoded, TSPY-like/SET/nucleosome assembly protein-1 superfamily. The encoded protein is localized to the nucleolus where it functions in chromatin remodeling and as an inhibitor of cell-cycle progression. This protein may play a role in the suppression of tumor growth.

**Antigen:** Synthetic peptide of human TSPYL2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitously expressed, with highest levels in brain, testis and heart, and lowest levels in liver and pancreas.