

**Anti-Human/Mouse CDC34 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09522**

---

3th Edition

**Description:**The protein encoded by this gene is a member of the ubiquitin-conjugating enzyme family. Ubiquitin-conjugating enzyme catalyzes the covalent attachment of ubiquitin to other proteins. This protein is a part of the large multiprotein complex, which is required for ubiquitin-mediated degradation of cell cycle G1 regulators, and for the initiation of DNA replication.

**Antigen:**Synthetic peptide of human CDC34

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in testes during spermatogenesis to regulate repression of cAMP-induced transcription.