

**Anti-Human/Mouse/Rat FRK Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09541**

---

3th Edition

**Description:**The protein encoded by this gene belongs to the TYR family of protein kinases. This tyrosine kinase is a nuclear protein and may function during G1 and S phase of the cell cycle and suppress growth.

**Antigen:**Recombinant protein of human FRK

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Restricted to cells lines derived from tissues of lymphoid, brain, breast, colon and bladder origin.