

Anti-Human/Mouse/Rat FSTL1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09544

3th Edition

Description: This gene encodes a protein with similarity to follistatin, an activin-binding protein. It contains an FS module, a follistatin-like sequence containing 10 conserved cysteine residues. This gene product is thought to be an autoantigen associated with rheumatoid arthritis.

Antigen: Recombinant protein of human FSTL1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Overexpressed in synovial tissues from rheumatoid arthritis.