

Anti-Human/Mouse CAMTA1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09551

3th Edition

Description: Calmodulin-binding transcription activator 1 is a protein that in humans is encoded by the CAMTA1 gene. Normally expressed in non-neoplastic adult central nervous system tissues: detected in whole brain, cerebellum, brain cortex, occipital lobe, frontal lobe, temporal lobe, putamen. Expression levels are low in oligodendroglial tumors, and are reduced by half in oligodendroglioma and astrocytoma cases with 1p loss of heterozygosity. Detected in neuroblastic-type cultured neuroblastoma cells.

Antigen: Synthetic peptide of human CAMTA1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Normally expressed in non-neoplastic adult central nervous system tissues: detected in whole brain, cerebellum, brain cortex, occipital lobe, frontal lobe, temporal lobe, putamen. Expression levels are low in oligodendroglial tumors, and are reduced by half in oligodendroglioma and astrocytoma cases with 1p loss of heterozygosity. Detected in neuroblastic-type cultured neuroblastoma cells. Expressed in heart and kidney.