

Anti-Human CARD17 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09607

3th Edition

Description: Regulator of procaspase-1/CASP1 activation implicated in the regulation of the proteolytic maturation of pro-IL-1 β /IL1B and its release during inflammation. Inhibits the release of IL1B in response to LPS in monocytes. However, unlike CASP1, do not induce NF-kappa-B activation.

Antigen: Synthetic peptide of human CARD17

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1 μ g/ml.

Specificity: Ubiquitous.